

join a **microscopy** revolution







revolution laser microscopy:

Andor Revolution - A Family of Products

Andor Revolution provides a framework for our range of laser-illuminated microscopy solutions, including spinning disk confocal, TIRFM and photo-bleach and activation modalities. Our spinning disk, live cell confocal instrument, illustrated below, combines Andor's award winning iXon EMCCD cameras, solid-state laser combiner and the renowned Yokogawa CSU confocal scanner "inside". Andor has established a global distribution agreement with Yokogawa Electric Corporation and harnessed our opto-electronics engineering and manufacturing knowledge to create powerful live cell confocal solutions. Andor's direct market presence and unique solutions bring you unrivalled performance, product understanding and support.



Revolution encompasses a range of components, both hardware and software that fit seamlessly together creating a complete laser microscopy solutions. An open flexible component focus also allows us to provide key pieces of hardware stand-alone. We want to tailor solutions matched to your needs. At the core of Revolution systems is Andor iQ, a multi-dimensional imaging software which synchronizes iXon EMCCD cameras with the CSU confocal spinning disk and other key hardware components such as Piezo Z100 z-stage and our Laser Combiner with AOTF for rapid laser line selection.

revolutionary technology:



Confocal Dual Spinning Disk

Confocal Dual Spinning Disk

Revolution's confocal dual spinning disk technology provides an ideal platform for high speed, high SNR imaging, with low bleach rate and low photoxicity.

- Isolated detection channel with extremely low background - optimal for EMCCD technology
- Efficient illumination coupling via a microlens array (~60%)
- Low peak power by parallel illumination of ~1000 scanning laser points
- Efficient parallel detection with highest QE 90% peak
- True broadband operation excitation 400-650nm - detection 420-750nm
- CSU optimization includes:
 - Apochromatic detector channel optics
 - Focusing single and dual camera adapters
 - Semrock dichroic filters optimized for selected laser lines

vailable Live Cell Confocal Upgrades

Andor's Revolution range includes our laser spinning disk live cell confocal microscopy systems utilizing the Yokogawa CSU. We have created a range of upgrade solutions for Perkin Elmer UltraView and Visitech system users to bring a new lease of life to legacy systems.

TIRF and Photo-bleaching

Andor Revolution technology is ideally suited to use in other modalities of fluorescence imaging, including TIRFM and photo-bleaching/activation. The Andor laser combiner (ALC) offers unique functionality for creation of multi-modal imaging systems.



Prism-type TIRFM (LHS) and objective-type TIRFM (RHS)

TIRFM is a well-established technique for live cell imaging making use of total internal reflection (TIR). The electromagnetic field of the TIR light (the evanescent field) extends into the sample beyond the interface by some one to two hundred nanometers into the medium of lower refractive index. Since the field decreases exponentially along the z-axis, only a very thin section of the specimen undergoes fluorescence excitation. TIRFM produces extremely high contrast, low background images, ideal for detection with our iXon EMCCD cameras and is a valuable tool in the study of e.g. cellular membrane traffic, surface receptor activity and cytoskeleton dynamics as well as in single molecule and molecular machine studies.

Andor Revolution components provide:

- Matched illuminator and detector for maximum performance
- Flexible fiber-coupled 400-650 nm excitation compatible with all TIRF platforms.
- Multiple fiber outputs for confocal, TIRF and Photobleach.
- Integral AOTF microsecond wavelength, shuttering and intensity control
- Active blanking precise synchronization and > 80 dB isolation.

microscopy products:

iXon[™]+ providing single photon sensitivity



Model	Features
• DU-680	128 x 128 pixels; 24 x 24µm; QE >90% for high speed applications e.g Ca²+ flux
• DU-885	1002 x 1004 pixels; 8 x 8µm; QE >60%, 35 MHz readout, high resolution intra- cellular imaging
• DU-897	512 x 512 pixels; 16 x 16µm; QE >90%, the ultimate combination of sensitivity and speed
• DU-888	1024 x 1024 pixels: 13 x 13um: QE >90%.

resolution and exceptional sensitivity combined



evolutio

2nio

All iXon detectors benefit from: a) minimized dark current from unparalleled –90°C TE cooling – especially critical for confocal imaging where experimental backgrounds are minimal, b) Industry fastest vertical shift rates (< 2µs/row) - Vital to avoid Scan Lines with both CSU 10 and CSU 22 operation

csu confocal spinning disk unit



Features

- 1500-5000 rpm disk, giving 300-1000 fps frame rate
- Easy synchronization of exposure to disk speed
- Manual and motorized filter changer models available
- FWHM Z section <=1.2µm
- Semrock dichroic and emission filter pairing provides optimum SNR
- Filter wheel mounting for fluorescence and DIC analyzer
- High performance custom Semrock dichroic mirrors; quad, triple and dual line combinations
- Custom optics fitted to camera adapters for higher throughput, chromatic aberration and background noise reduction

andor laser combiner solid state & gas laser lines,



Solid State Range

- 405 50 mW
- 445 40 mW
- 488 20, 25, 50 mW
- 513 20, 50 mW
- 532 20, 50, 75 mW
- 561 25, 50, 75 mW
- 635 25 mW

Note: 400 series model allows up to 4 solid state lasers (5 lines with a gas laser), while the soon to be announced 600 series will support up to 6 solid state lasers

Features

- Compact 4U 19" rack mount enclosure
- Quiet, cool, low power consumption
- Long life solid state lasers with excellent stability (typical +/- 0.4% pk to pk)
- Integral AOTF for laser selection. intensity
 and blanking control
- Excellent isolation and cross-talk (>80dB)
- ALC allows Gas laser integration e.g. Ar ion or mixed gas

precision controller unit flexible expansion



PCU100

- Provides system synchronization
- Houses AOTF driver for ALC coupling
- 16 bit Analogue output for direct piezo Z control
- 8 digital inputs for sensing external events
- 8 digital outputs for triggering and expansion
- Enables Revolution to slave or master.
- All under iQ software control

PCU NEOS

- Up to 5x 16 bit analogue outputs for external AOTF control
- NEOS or AAEO AOTF support
- Optional active blanking driver for AOTF shuttering
- 8 digital inputs for sensing events
- 8 digital outputs for triggering and expansion
- Information available for third party software support

motorized xy & z control piezo z / motorized xy



XY Stage	Z Stage	Z Objective	
 Open and closed loop stages for your preferred 	 Piezo Z stage with 100, 200 and 400µm ranges 	 PI PIFOC P721-100µm or P725-400µm travel 	
microscope	 1nm repeatability 	• Up to 20 Z sections per	
• Typical travel >100x75mm,	 2-8 nm resolution 	second	
with 0.02µm resolution	(depends on range)	 1.25nm resolution 	
 30mm/sec travel speed 	Up to 100 Z sections per	Analogue or digital control	
Perform multi-field scans	second	Oil and water objectives	
for 6D imaging	 Slide holder 3"x1" 		

Create 4D mosaics using iQ software

support of Sutter

Multi-wheel Sutter controllers available

Smartshutter

• 35mm petri dish holder

filter wheels & splitters emission discrimination



TR-FW10	TR-FW06	Field Splitte
 10 position, ø25 filterwheel with 40ms gwitching time 	• 6 position, ø25 filterwheel,with 100 ms	Optosplit II us filters
 Matched with CSU or microscope "C"coupler 	 Matched with CSU or microscope "C" coupler 	 Field splitting tail simultaneous imaging
 Infinity optical path for optimum filter 	 Infinity optical path for optimum filter 	 Matched to C microscope "
performanceExternal controller allows	performanceMAC5000 external	 DualView and also supporte

controller expandable for other Ludl components e.g. stages

- es 25 mm
- for 2 color
- SU or C" coupler
- QuadView iso supported

software optimized imaging software & choice



iQ

Andor iQ is our flagship live cell imaging software, designed with flexibility and power in mind. iQ - image and quantify - provides optimized control of Andor's award winning iXon^{EM}+ EMCCD cameras and other hardware for bio-imaging.

Imaris v7.0 and VisBio Bio-formats import Andor Multi-TIFF image files.

MatLab - We are a member of Mathworks connections program and provide Revolution data access tools for MatLab users.

Micro-Manager and Image-J -Revolution components will soon be compatible with the Micro-Manager open source device adapter standard.

iQ Workstation

- Dell Precision workstation
- Preconfigured with iQ and Control hardware
- RAID-0 disk configuration
- High performance graphics with OpenGL
- Guaranteed compatibility and maximum frame rate performance

camera port adapters single and dual port



Single Port Coupler

- Magnification X = 1.0, 1.2, 1.5, 2.0 other by request
- Compatible with filter wheel or CSU+filter wheel

Andor couplers are designed to ensure optimum throughput, minimal aberrations, magnification options and flexible detection configurations. Our couplers match our CSU enhancements and can also be configured with a broad range of imaging optics including microscope C-mount, 'C' and 'CS' mount lenses.

Dual Port Coupler



Jolut

- Magnification X = 1.0, 1.2, 1.5, 2.0 other by request
- Compatible with 'C'/'CS' lenses, C-mount, filter wheel or CSU+filter wheel
- Magnification select in each arm
- Mounting cassettes 'store adjustment' simplifying exchange and minimizing readjustment

microscopes







Diverse Compatibility

Revolution products are compatible with modern infinity corrected microscopes from Leica, Nikon, Olympus or Zeiss to meet your preferred configuration. NB For laser safety a 100% side port switch is required. We have recently introduced support for TILL's iMIC - modular automated microscope system.

Microscopes need environmental control to maintain living conditions for mammalian cell preparations. We recommend temperature and CO² control incubators from LIS (Switzerland). Thermal and mechanical factors in the microscope mechanics or environment can impact the observation of living specimens when image focus is affected by drift or vibration. Therefore we recommend that systems are configured on vibration isolation tables. An important recent innovation in microscope design is focus drift correction (FDC). FDC corrects focus drift, which can severely affect time-lapse movies and is especially important in confocal and TIRF imaging, where focus and illumination are tightly constrained. Nikon and Olympus have introduced FDC solutions, while ASI (Oregon, USA) offer a product they call "CRIFF" which is adaptable to most microscopes. Andor iQ supports the operation of these devices, to provide long-term drift free imaging from Revolution live cell imaging solutions.

microscopy applications:

Revolution components and systems benefit from common features across applications. Whether in confocal, photo-bleach/activation or TIRFM imaging this means state of the art performance, flexible component-oriented supply and an open-minded approach to systems and technically demanding solutions. As we continue to provide wide-ranging support for key hardware components, we also try to increase your choice through expanding our software compatibility. We are actively involved in making the bridge from iQ to MatLab, Imaris or Image-J and Micro-Manager as smooth as possible.

Applications include:

- Calcium Signaling
- Multi-channel GFP Imaging
- Developmental and cell cycle studies
- Fluorescence translocation and Colocalization
- FRET imaging
- Membrane trafficking
- Cytoskeleton dynamics
- Single molecule studies
- Molecular machines

System Features

- Ultimate sensitivity from combination of iXon[™]+ EMCCD and TIRF illumination or CSU spinning disk.
- Fast frame rates (up to 1000 sec⁻¹) with perfect synchronization.
- Reduced photo bleaching particularly vital for nonratiometric dyes.
- Reduced phototoxicity of living cells.
- Study live specimens with reduced fluorophore concentrations or expression levels.
- Minimal perturbation of physiological events.
- Synchronization of electrophysiology and photolysis with the Precision Control unit.
- Rapid optical sectioning with Piezo Z100 for 4D imaging up to 100 sections sec⁻¹.
- Detection limit down to single GFP molecules in parallel at supra-video frame rates.
- Powerful visualization, image analysis and tracking tools in iQ software.
- Detect and track many individual vesicles in parallel.
- Lower detection limits down to single co-localized pairs, at supra-video frame rates.
- Align and merge optically split images in iQ.
- Powerful and comprehensive FRET analysis and tracking.



Courtesy Dr lain Johnson, Invitrogen, Eugene, OR. Cultured neurons were loaded with DiSBAC2(3), a voltage sensitive dye and visualized with the Revolution 488. Exposure time of the IXon 887-BV was 100ms and EM Gain 200. The figure shows a film strip from Andor iQ of maximum intensity projection time series. Time point 4 is taken almost simultaneously with the addition of KCl, which depolarizes the cells and shows a near instantaneous rise in DISBAC2(3) signal. ΔF of almost 45% is observed in this example. Imaged at the 3D Microscopy of Imaging Cells Course, Vancouver 2005.



Courtesy Dr James E. Bear, Lineberger Comprehensive Cancer Center, NC. Rat2 fibroblast expressing Pod-GFP (green) stained with phalloidin (red).













"With iXon EMCCD cameras we are seeing a revolution, the revolution is seeing (single molecules in action). The synergy of iXon and Yokogawa CSU technology will open new exciting possibilities in cellular imaging"

Dr Derek K Toomre Assistant Professor of Cell Biology Member of Ludwig Institute for Cancer Research Yale University School of Medicine





Embryo Cell Division using the Andor Revolution XD system. Courtesy of Paula Sampaio, University of Porto Rua Campo Alegre, Portugal

Andor **overview:**

Andor Technology strives to bring innovative solutions to the demanding needs of researchers and manufacturers. We are a reliable source for the researcher in the lab as well as a trusted supplier to more than half of the top twenty analytical instrumentation companies in the world. We understand their requirements and share an in-depth knowledge of their applications. Whether you require single or multiple components, off-the-shelf or a custom-designed solution, Andor is the logical choice, enabling greater efficiencies, lower costs and better results.



Andor Headquarters, Belfast



www.andor.com

Andor Technology Worldwide: +44 28 9023 7126 Andor Technology USA: (860) 290-9211 Andor Technology Japan: 81-3-3511 0659 Andor Technology China: +86-10-5129-4977



join the **revolution**





March 2007